


Dental Radiology range - X-Mind



Length of wall bracket (meters)	Arm extension (meters)
0.40	1.43
0.80	1.84
1.10	2.12



X-Mind AC
 Generator: low frequency (50Hz)
 Minimum exposure time: 0.08 seconds
 Voltage: 70KV
 Intensity: 8mA

X-Mind DC
 High-frequency generator (200Hz): constant voltage
 Minimum exposure time: 0.02 seconds
 Reduction of exposure time: 35% compared with X-Mind AC
 Voltage: 70KV or 60KV
 Intensity: 4mA or 8mA
 Reduction of radiation dose: 80% with a digital sensor

X-Mind
mobile

X-Mind
wall mounted